


PRODUCT ID	E-1178	<p align="center">Structure</p> 
PRODUCT	Androsin/草夹竹桃苷	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	15090122	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2015-09-01	
DATE OF REPORT	2015-09-01	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₁₅H₂₀O₈

MOLECULAR WEIGHT(MW)

328.317

CAS NUMBER

531-28-2

CHEMICAL FAMILY

FROM

Figwortflower Picrorhiza Rhizome

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK49—2006	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK49—2006	NA	ND
Water	Q/TFFK49—2006	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Androsin **Chrom-Column:** C18
Lot Number: 15090122 **Mobile Phase:** ACN:0.1%H3PO4=10:90
Detector: UV **Wavelength:** 268nm
Detect Time: 2015-09-01 **Inject Vol:** 20ul
Analyst: Shuzhilin **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	9.128	1641.174	22505.500	0.4477
2	13.355	26003.219	5004673.500	99.5523