



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0926	Structure 
PRODUCT	Salvianolic acid D/丹酚酸D	
STANDARD TYPE	Qualitative Reference Material/ 含量测定	
LOT NUMBER	17032822	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2017-03-28	
DATE OF REPORT	2017-03-28	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₂₀H₁₈O₁₀

MOLECULAR WEIGHT(MW)

418.35

CAS NUMBER

142998-47-8

CHEMICAL FAMILY

Phenylpropanoids 苯丙素类

FROM

Salvia miltiorrhiza Bge. 丹参

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK170—2007	NLT ≥96.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK170—2007	NA	ND
Water	Q/TFFK170—2007	≤2%	Conforms
Appearance		Pale Yellow Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Salviannolic acid D **Chrom-Column:** C18
Lot Number: 17032822 **Mobile Phase:** ACN:0.1%H₃PO₄=26:74
Detector: UV **Wavelength:** 320nm
Detect Time: 2017-03-28 **Inject Vol:** 20ul
Analyst: shuzhilin **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	6.535	356.104	4622.720	0.0377
2	7.044	892084.938	11882830.000	96.7983
3	15.539	7720.992	207654.469	1.6916
4	17.090	5743.675	180761.641	1.4725