



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0249	Structure 
PRODUCT	Pseudoginsenoside F11 /拟人参皂苷F11	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	18012922	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-01-29	
DATE OF REPORT	2018-01-29	

OTHER NAME	NA
CHEMICAL FORMULA	C42H72O14
MOLECULAR WEIGHT(MW)	801.01268
CAS NUMBER	69884-00-0
CHEMICAL FAMILY	Triterpenes

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK247—2007	NLT ≥98.0%	100.0%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK247—2007	NA	ND
Water	Q/TFFK247—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst

Checker

Manager

Product Analysis Report

Product Name: Pseudoginsenoside F11 Chrom-Column: C18

Lot Number: 18012922

Mobile Phase: ACN:H2O (33: 67)

Detect Time: 2018-1-29

Detector: ELSD

Inject Vol: 20 μ L


Tube Temp: 105 $^{\circ}$ C

Flow: 1.0ml/min

Gas Flow: 2.5 L/min

Sample MeOH

Impactor: Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	6.665	720078.375	7705396.000	100.0000