



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-1143	Structure 
PRODUCT	Vincetoxicoside B/田基黄苷	
STANDARD TYPE	Quantitative Reference Standard/ 含量测定	
LOT NUMBER	15101222	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2015-10-12	
DATE OF REPORT	2015-10-12	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA C₂₁H₂₀O₁₁

MOLECULAR WEIGHT(MW) 448.38

CAS NUMBER 22007-72-3

CHEMICAL FAMILY

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK21—2004	NLT ≥97.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK21—2004	NA	ND
Water	Q/TFFK21—2004	≤2%	Conforms
Appearance		Pale Yellow Powder	Conforms

STORAGE CONDITIONS

STORAGE Refer to the label
SHELF LIFE 2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Vincetoxicoside B **Chrom-Column:** C18
Lot Number: 15101222 **Mobile Phase:** MeOH:0.1%H3PO4(45:55)
Detector: UV **Wavelength:** 255nm
Detect Time: 2015-10-12 **Inject Vol:** 20ul
Analyst: zhufangfang **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	18.090	402.532	21418.600	0.1807
2	29.373	141123.469	11504874.000	97.0768
3	35.507	2766.967	218862.422	1.8467
4	37.507	673.674	55806.590	0.4709
5	42.107	72.535	7131.646	0.0602
6	49.573	318.178	35947.301	0.3033
7	56.107	70.500	7266.400	0.0613