



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-1176	<p align="center">Structure</p> 
PRODUCT	2-Acetylacteoside /2-乙酰基毛蕊花糖苷	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	15070832	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2015-07-08	
DATE OF REPORT	2015-07-08	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₃₁H₃₈O₁₆

MOLECULAR WEIGHT(MV)

666.624

CAS NUMBER

94492-24-7

CHEMICAL FAMILY

FROM

Cistanche salsa

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK1—2003	NLT ≥97.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK1—2003	NA	ND
Water	Q/TFFK1—2003	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: 2-Acetylacteoside **Chrom-Column:** C18
Lot Number: 15070832 **Mobile Phase:** MeOH:0.1%H₃PO₄=4
Detector: UV **Wavelength:** 330nm
Detect Time: 2015-07-08 **Inject Vol:** 20ul
Analyst: shuzhilin **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	4.842	1185.824	14924.650	0.1539
2	5.701	216.093	4693.800	0.0484
3	6.485	450.303	7934.313	0.0818
4	7.252	494182.813	9456593.000	97.4901
5	8.019	794.297	16710.377	0.1723
6	10.103	161.410	3237.590	0.0334
7	10.670	742.841	15832.822	0.1632
8	11.237	2037.273	85715.992	0.8837
9	14.506	228.750	13924.700	0.1436
10	18.908	1176.385	80484.797	0.8297