



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-2319	<p align="center">Structure</p> 
PRODUCT	7,2'-Dihydroxy-3',4'-dimethoxyisoflavan-7-O-beta-D-glucopyranoside /7,2' -二羟基-3',4' -二甲氧基-异黄烷-7-O-β-D-葡萄糖苷	
STANDARD TYPE	Quantitative Reference Standard/ 含量测定	
LOT NUMBER	15071522	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2015-07-15	
DATE OF REPORT	2015-07-15	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₂₃H₂₈O₁₀

MOLECULAR WEIGHT(MW)

464.46

CAS NUMBER

94367-43-8

CHEMICAL FAMILY

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK417—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK417—2007	NA	ND
Water	Q/TFFK417—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst


Checker

Manager

Product Analysis Report

Product Name: 7,2-Dimethoxy-3',4'-dimethoxyisoflavan-7-O-β-D-glucopyranoside
Lot Number: 15071522
Temperature: 25°C
Detect Time: 2015-07-15
Analyst: shuzhilin

Chrom-Column: C18
Mobile Phase: MeOH:H2O (grads)
Wavelength: 280nm
Inject Vol: 20μL
Flow: 1ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	11.650	496.004	7326.549	0.1532
2	14.617	341390.219	4765097.000	99.6616
3	21.582	658.625	8852.551	0.1852

	grads	
Time/min	0	30
MeOH	40	80
H2O	60	20