


PRODUCT ID	E-1117	<p style="text-align: center;">Structure</p> 
PRODUCT	Ethyl 4'-hydroxy-3'-methoxycinnamate/阿魏酸乙酯	
STANDARD TYPE	Quantitative Reference Standard/含量测定	
LOT NUMBER	15060522	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2015-06-05	
DATE OF REPORT	2015-06-05	

CHEMICAL NAME
OTHER NAME
CHEMICAL FORMULA

 C₁₂H₁₄O₄
MOLECULAR WEIGHT(MW)

222.24

CAS NUMBER

4046-02-0

CHEMICAL FAMILY
FROM
ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK42—2006	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK42—2006	NA	ND
Water	Q/TFFK42—2006	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

 STORAGE
 SHELF LIFE

 Refer to the label
 2 Years

Analyst


Checker

Manager

Product Analysis Report

Product Name: Ethyl 4'-hydroxy-3'-methoxycinnamate
Lot Number: 15060522
Detector: UV
Detect Time: 2015-06-05
Analyst: shuzhilin

Chrom-Column: C18
Mobile Phase: MeOH:H2O=60:40
Wavelength: 320nm
Inject Vol: 20ul
Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	9.870	349339.313	8071421.500	100.0000