



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

MSN:Tauto_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-1255	Structure 
PRODUCT	Piceatannol/白皮杉醇	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	15040823	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2015-04-08	
DATE OF REPORT	2015-04-08	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₁₄H₁₂O₄

MOLECULAR WEIGHT(MW)

244.25

CAS NUMBER

10083-24-6

CHEMICAL FAMILY

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK248—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK248—2007	NA	ND
Water	Q/TFFK248—2007	≤2%	Conforms
Appearance		Pale Yellow Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Piceatannol **Chrom-Column:** C18
Lot Number: 15040823 **Mobile Phase:** ACN:0.1% H3PO4(20:80)
Detector: UV **Wavelength:** 319nm
Detect Time: 2015-04-08 **Inject Vol:** 20ul
Analyst: shuzhilin **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	17.891	240368.422	6636065.500	99.7645
2	32.264	184.846	8301.969	0.1248
3	34.732	156.585	7365.100	0.1107

0)

