


PRODUCT ID	E-1058	<p>Structure</p> 
PRODUCT	3,29-Dibenzoyl rarounitriol /3,29-二苯甲酰基栝楼仁三醇	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	14090325	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2014-09-03	
DATE OF REPORT	2014-09-03	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C44H58O5

MOLECULAR WEIGHT(MW)

666.93

CAS NUMBER

873001-54-8

CHEMICAL FAMILY

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK762—2010	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK762—2010	NA	ND
Water	Q/TFFK762—2010	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: 3,29-Dibenzoyl rarounitriol Chrom-Column: C18
Lot Number: 14090325 Mobile Phase: MeOH:H2O=93:7
Temperature: 25°C Wavelength: 230nm
Detect Time: 2014-09-03 Inject Vol: 20.0 µL
Analyst: shuzhilin Flow: 1.00 ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	31.523	456777.750	30229840.000	98.6563
2	34.643	3497.252	411720.563	1.3437