



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

MSN:Tauto_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-2468	Structure 
PRODUCT	Ginsenoside F3/人参皂苷F3	
STANDARD TYPE	Qualitative Reference Material/ 含量测定	
LOT NUMBER	15082024	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2015-08-20	
DATE OF REPORT	2015-08-20	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₄₁H₇₀O₁₃

MOLECULAR WEIGHT(MW)

770.99

CAS NUMBER

62025-50-7

CHEMICAL FAMILY

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK122—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK122—2007	NA	ND
Water	Q/TFFK122—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst


Checker

Manager

Product Analysis Report

Product Name: Ginsenoside F3
Lot Number: 15082024
Temperature: 25°C
Detect Time: 2015-08-20
Analyst: shuzhilin

Chrom-Column: C18
Mobile Phase: ACN:H2O(35:65)
Wavelength: 203nm
Inject Vol: 20ul
Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	7.778	821.898	15691.727	0.3446
2	8.527	308354.906	4538298.500	99.6554