


PRODUCT ID	E-0835	<p>Structure</p> 
PRODUCT	Liriope muscari baily saponins C/短 葶山麦冬皂苷C	
STANDARD TYPE	Quantitative Reference Standard/ 含 量测定	
LOT NUMBER	16022522	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2016-02-25	
DATE OF REPORT	2016-02-25	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₄₄H₇₀O₁₆

MOLECULAR WEIGHT(MV)

855.02

CAS NUMBER

87480-46-4

CHEMICAL FAMILY

Steroids甾体类

FROM

麦冬*Ophiopogon japonicus*

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK320—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK320—2007	NA	ND
Water	Q/TFFK320—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Liriope muscari baily saponins C Chrom-Column: Arcus C18
Lot Number: 16022522 Mobile Phase: ACN:H2O(46: 54)
Detect Time: 2016-2-25 Detector: ELSD
Inject Vol: 20 μ L Tube Temp: 105 $^{\circ}$ C
Flow: 1.0ml/min Gas Flow: 2.5L/min
Analyst: zhufangfang Impactor: Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	8.473	149125.500	2000098.875	100.0000