



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

MSN:Tauto_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0901	<p>Structure</p>  <p>86831-54-1 Rhodiosin 红景天素</p>
PRODUCT	Rhodiosin/红景天素	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	16011931	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2016-01-19	
DATE OF REPORT	2016-01-19	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₂₇H₃₀O₁₆

MOLECULAR WEIGHT(MW)

610.52

CAS NUMBER

86831-54-1

CHEMICAL FAMILY

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK165—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK165—2007	NA	ND
Water	Q/TFFK165—2007	≤2%	Conforms
Appearance		plae yellow powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst


Checker

Manager

Product Analysis Report

Product Name: Rhodiosin
Lot Number: 16011931
Temperature: 25°C
Detect Time: 2016-01-19
Analyst: shuzhilin

Chrom-Column: C18
Mobile Phase: ACN:0.1%H3pPO4=26:74
Wavelength: 325nm
Inject Vol: 20ul
Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	7.973	593.735	11403.300	0.2931
2	10.598	177895.984	3879675.500	99.7069