



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-2574	<p align="center">Structure</p> 
PRODUCT	Helenien/叶黄素棕榈酸酯(叶黄素酯)	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	16101332	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2016-10-13	
DATE OF REPORT	2016-10-13	

CHEMICAL NAME

OTHER NAME

叶黄素二棕榈酸酯、堆心菊素

CHEMICAL FORMULA

C72H116O4

MOLECULAR WEIGHT(MW)

1045.69

CAS NUMBER

547-17-1

CHEMICAL FAMILY

aliphatic ester类胡萝卜素脂肪酸酯

FROM

万寿菊Tagetes erecta L

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK680—2006	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK680—2006	NA	ND
Water	Q/TFFK680—2006	≤2%	Conforms
Appearance		Yellow Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst

Checker

Manager

