



同田生物  
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

<b>PRODUCT ID</b>	E-3161	<b>Structure</b> 
<b>PRODUCT</b>	D-Mannitol/D-甘露醇	
<b>STANDARD TYPE</b>	Quantitative Reference Material (含量测定)	
<b>LOT NUMBER</b>	17111324	
<b>ASSAY METHOD</b>	HPLC,NMR	
<b>DATE OF SAMPLE</b>	2017-11-13	
<b>DATE OF REPORT</b>	2017-11-13	

**OTHER NAME**

**CHEMICAL FORMULA**

C<sub>6</sub>H<sub>14</sub>O<sub>6</sub>

**MOLECULAR WEIGHT(MW)**

182.17

**CAS NUMBER**

69-65-8

**CHEMICAL FAMILY**

Aliphatic

**FROM**

**ANALYTICAL RESULT**

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK747—2010	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK747—2010	NA	ND
Water	Q/TFFK747—2010	≤2%	Conforms
Appearance		white Powder	Conforms

**STORAGE CONDITIONS**

STORAGE  
SHELF LIFE

Refer to the label  
2 Years


Analyst

Checker

Manager

# Product Analysis Report

**Product Name:** D-Mannitol      **Chrom-Column:** Arcus C18  
**Lot Number:** 17111324      **Mobile Phase:** H2O  
**Detect Time:** 13-Nov-17      **Detector:** ELSD  
**Inject Vol:** 20 $\mu$ L      **Tube Temp:** 105 $^{\circ}$ C  
**Flow:** 1.0ml/min      **Gas Flow:** 2.5L/min  
**Analyst:** zhanghong      **Impactor:** Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	3.265	612682.688	4017500.750	100.0000