



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com

export2@tautobiotech.com

PRODUCT ID	E-0688	<p align="center">Structure</p> 
PRODUCT	Akebiasaponin D/川续断皂苷 VI/ 木通皂苷 D	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17071441	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2017-07-14	
DATE OF REPORT	2017-07-14	

OTHER NAME Asperosaponin VI

CHEMICAL FORMULA C47H76O18

MOLECULAR WEIGHT(MW) 929.1

CAS NUMBER 39524-08-8

CHEMICAL FAMILY

FROM 木通Akebia quinata (Thunb.) Decne

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK539—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK539—2007	NA	ND
Water	Q/TFFK539—2007	≤2%	Conforms
Appearance		white Powder	Conforms

STORAGE CONDITIONS

STORAGE Refer to the label

SHELF LIFE 2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Akebiasaponin D **Chrom-Column:** C18
Lot Number: 17071441 **Mobile Phase:** ACN:H2O (30:70)
Detect Time: 2017-7-14 **Detector:** ELSD
Inject Vol: 20.0 μ L **Tube Temp:** 105 $^{\circ}$ C
Flow: 1.0ml/min **Gas Flow:** 2.5 L/min
Analyst: zhanghong **Impactor:** Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	15.432	148153.828	3182416.250	100.0000