



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

MSN:Tauto_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-3069	<p style="text-align: center;">Structure</p> 
PRODUCT	Astragaloside VI/黄芪皂苷VI	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17072021	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2017-07-20	
DATE OF REPORT	2017-07-20	

OTHER NAME

CHEMICAL FORMULA C47H78O19

MOLECULAR WEIGHT(MW) 946

CAS NUMBER 84687-45-6

CHEMICAL FAMILY Triterpenoid三萜类

FROM 膜荚黄芪Astragalus membranaceus (Fisch.) Bunge

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK143—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK143—2007	NA	ND
Water	Q/TFFK143—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE Refer to the label

SHELF LIFE 2 Years

Analyst

Checker


Manager

Product Analysis Report

Product Name: Astragaloside VI Chrom-Column: C18

Lot Number: 17072021 Mobile Phase: ACN:H2O=32: 68

Detect Time: 2017-7-20 Detector: ELSD
Inject Vol: 20 μ L Tube Temp: 105 $^{\circ}$ C
Flow: 1.0ml/min Gas Flow: 2.5 L/min
Analyst: zhanghong Impactor: Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	16.265	264872.469	4716500.000	100.0000