



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0689	Structure 
PRODUCT	Corilagin/柯里拉京	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17092723	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2017-09-27	
DATE OF REPORT	2017-09-27	

OTHER NAME

CHEMICAL FORMULA

C₂₇H₂₂O₁₈

MOLECULAR WEIGHT(MW)

634.4528

CAS NUMBER

23094-69-1

CHEMICAL FAMILY

Phenols酚类

FROM

大戟科油柑属植物叶下珠Phyllanthus urinaria L. 全草

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK565—2008	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK565—2008	NA	ND
Water	Q/TFFK565—2008	≤2%	Conforms
Appearance		White powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Corilagin Chrom-Column: C18
Lot Number: 17092723 Mobile Phase: ACN:0.1%H₃PO₄=15:85
Detector: UV Wavelength: 278nm
Detect Time: 2017-9-27 Inject Vol: 20ul
Analyst: shuzhilin Flow: 1.0ml/min



NO.	Ret. Time min	Height mAU	Area mAU*min	Rel. Area%
1	7.577	463095.594	5489941.500	100.0000