



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

MSN:Tauto_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-3089	<p align="center">Structure</p> 
PRODUCT	Melittoside/蜜利特昔	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17071023	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2017-07-10	
DATE OF REPORT	2017-07-10	

OTHER NAME

CHEMICAL FORMULA

C₂₁H₃₂O₁₅

MOLECULAR WEIGHT(MW)

524.472

CAS NUMBER

19467-03-9

CHEMICAL FAMILY

Iridoids环烯醚萜类

FROM

Rehmannia glutinosa (Gaetn.) Libosch. ex Fisch. et Mey.
地黄

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK231—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK231—2007	NA	ND
Water	Q/TFFK231—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Melittoside Chrom-Column: C18
Lot Number: 17071023 Mobile Phase: ACN:0.1%H3PO4=5: 95
Detector: UV Wavelength: 205nm
Detect Time: 2017-07-10 Inject Vol: 20ul
Analyst: shuzhilin Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	4.184	1175.323	5773.900	0.0297
2	8.786	195.471	3241.200	0.0167
3	9.887	783.068	13365.993	0.0687
4	11.237	990907.313	19422518.000	99.8849