



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-2567	Structure 
PRODUCT	Uracil/尿嘧啶	
STANDARD TYPE	Quantitative Reference Standard/ 含量测定	
LOT NUMBER	17051133	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2017-05-11	
DATE OF REPORT	2017-05-11	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₄H₄N₂O₂

MOLECULAR WEIGHT(MW)

112.09

CAS NUMBER

66-22-8

CHEMICAL FAMILY

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK159—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK159—2007	NA	ND
Water	Q/TFFK159—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Uracil **Chrom-Column:** Arcus C18
Lot Number: 17051133 **Mobile Phase:** 100% 0.02MKH2PO4(pH=3.0)
Temperature: 25°C **Wavelength:** 254nm
Detect Time: 2017-05-11 **Inject Vol:** 20μL
Analyst: zhufangfang **Flow:** 1ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	7.415	613675.688	10179926.000	100.0000