


Product Analysis Report

Product Name: Cinnamyl alcohol **Chrom-Column:** Arcus C18
Lot Number: 17092722 **Mobile Phase:** MeOH:0.1%H3PO4=52: 48
Temperature: 25°C **Wavelength:** 254nm
Detect Time: 2017-09-27 **Inject Vol:** 20μL
Analyst: zhanghong **Flow:** 1ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	7.373	424.489	6118.400	0.0375
2	10.490	849989.063	16289583.000	99.8891
3	15.098	79.902	1744.100	0.0107
4	22.532	276.078	10223.200	0.0627