



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-1198	Structure 
PRODUCT	Spicatoside B/山麦冬皂苷B	
STANDARD TYPE	Quantitative Reference Standard/ 含量测定	
LOT NUMBER	17040931	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2017-04-09	
DATE OF REPORT	2017-04-09	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₃₉H₆₂O₁₂

MOLECULAR WEIGHT(MW)

722.4

CAS NUMBER

87425-34-1

CHEMICAL FAMILY

Steroids甾体类

FROM

百合科植物山麦冬的干燥块根

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK133—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK133—2007	NA	ND
Water	Q/TFFK133—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Spicatoside B **Chrom-Column:** Arcus C18
Lot Number: 17040931 **Mobile Phase:** ACN:H2O (60:40)
Detect Time: 2017-4-9 **Detector:** ELSD
Inject Vol: 20μL **Tube Temp:** 105°C
Flow: 1.0ml/min **Gas Flow:** 3.0L/min
Analyst: zhufangfang **Impactor:** Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	8.998	368709.781	5421648.500	100.0000