



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-1155	Structure 
PRODUCT	3-dehydrotrametenolic acid/ 松苓新酸	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17081232	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2017-08-12	
DATE OF REPORT	2017-08-12	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₃₀H₄₆O₃

MOLECULAR WEIGHT(MW)

454.69

CAS NUMBER

29220-16-4

CHEMICAL FAMILY

FROM

茯苓 (Poria cocos)

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK724-2010	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK724-2010	NA	ND
Water	Q/TFFK724-2010	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: 3-dehydrotrametenolic aci **Chrom-Column:** C18
Lot Number: 17081232 **Mobile Phase:** ACN:0.1%H3PO4=80:20
Detector: UV **Wavelength:** 240nm
Detect Time: 2017-08-12 **Inject Vol:** 20ul
Analyst: shuzhilin **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	8.102	597.015	5974.100	0.0371
2	25.244	613636.625	16095200.000	99.9629