



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-3426	Structure 
PRODUCT	Stevioside/甜菊醇	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17111521	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2017-11-15	
DATE OF REPORT	2017-11-15	

OTHER NAME 甜叶菊苷元

CHEMICAL FORMULA C₂₀H₃₀O₃

MOLECULAR WEIGHT(MW) 318.45

CAS NUMBER 471-80-7

CHEMICAL FAMILY Diterpenoids

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK747—2010	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK747—2010	NA	ND
Water	Q/TFFK747—2010	≤2%	Conforms
Appearance		white Powder	Conforms

STORAGE CONDITIONS

STORAGE Refer to the label

SHELF LIFE 2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Stevioside **Chrom-Column:** Arcus C18
Lot Number: 17111521 **Mobile Phase:** ACN:H2O (40%-70%)
Detect Time: 15-Nov-17 **Detector:** ELSD
Inject Vol: 20 μ L **Tube Temp:** 110 $^{\circ}$ C
Flow: 1.0ml/min **Gas Flow:** 2.5L/min
Analyst: zhanghong **Impactor:** Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	21.932	119313.531	3147026.500	100.0000