



**同田生物**  
**TAUTO BIOTECH**

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

<b>PRODUCT ID</b>	E-1229	<p align="center"><b>Structure</b></p> 
<b>PRODUCT</b>	Didymine /香风草苷	
<b>STANDARD TYPE</b>	Qualitative Reference Material/含量测定	
<b>LOT NUMBER</b>	17072722	
<b>ASSAY METHOD</b>	HPLC,MS,NMR	
<b>DATE OF SAMPLE</b>	2017-07-27	
<b>DATE OF REPORT</b>	2017-07-27	

<b>OTHER NAME</b>	香风草甙
<b>CHEMICAL FORMULA</b>	C <sub>28</sub> H <sub>34</sub> O <sub>14</sub>
<b>MOLECULAR WEIGHT(MV)</b>	594.56
<b>CAS NUMBER</b>	14259-47-3
<b>CHEMICAL FAMILY</b>	Flavonoids黄酮类
<b>FROM</b>	断血流HERBA CLINOPODII

**ANALYTICAL RESULT**

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK1—2003	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK1—2003	NA	ND
Water	Q/TFFK1—2003	≤2%	Conforms
Appearance		White Powder	Conforms

**STORAGE CONDITIONS**

STORAGE	Refer to the label
SHELF LIFE	2 Years


**Analyst**

**Checker**

**Manager**

# Product Analysis Report

**Product Name:** Didymin      **Chrom-Column:** C18  
**Lot Number:** 17072722      **Mobile Phase:** ACN:0.1%H3PO4=28:72  
**Detector:** UV      **Wavelength:** 276nm  
**Detect Time:** 2017-07-27      **Inject Vol:** 20ul  
**Analyst:** shuzhilin      **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	12.296	3012.045	68486.500	0.5897
2	14.155	437205.438	11539781.000	99.3685
3	16.640	153.723	4855.672	0.0418