



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0828	<p align="center">Structure</p> 
PRODUCT	Narcissoside/异鼠李素-3-芸香糖甙/水仙苷	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17051121	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2017-05-11	
DATE OF REPORT	2017-05-11	

CHEMICAL NAME

OTHER NAME

Isorhamnetin 3-O-rutinoside

CHEMICAL FORMULA

C₂₈H₃₂O₁₆

MOLECULAR WEIGHT(M

624.54

CAS NUMBER

604-80-8

CHEMICAL FAMILY

Flavonoids黄酮类

FROM

Narcissus tazetta水仙

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK692—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK692—2007	NA	ND
Water	Q/TFFK692—2007	≤2%	Conforms
Appearance		yellow Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Narcissoside Chrom-Column: C18
Lot Number: 17051121 Mobile Phase: MeOH:ACN:0.1%H₃PO₄=15:15:70
Temperature: 25°C Wavelength: 260nm
Detect Time: 2017-05-11 Inject Vol: 20ul
Analyst: Shuzhilin Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	19.483	1242.000	31185.400	0.6629
2	21.476	166769.000	4665339.500	99.1634
3	23.277	90.000	3496.000	0.0743
4	24.444	110.000	4676.800	0.0994