



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-1033	Structure 
PRODUCT	6''-Feruloylspinosin/6''阿魏酰基斯皮诺素	
STANDARD TYPE	Quantitative Reference Standard/含量测定	
LOT NUMBER	18061421	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2018-06-14	
DATE OF REPORT	2018-06-14	

OTHER NAME	Spinosin B
CHEMICAL FORMULA	C ₃₈ H ₄₀ O ₁₈
MOLECULAR WEIGHT(MW)	784.71
CAS NUMBER	77690-92-7
CHEMICAL FAMILY	Flavonoids

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK774—2010	NLT ≥98.0%	99.8%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK774—2010	NA	ND
Water	Q/TFFK774—2010	≤2%	Conforms
Appearance		yellow Powder	Conforms

STORAGE CONDITIONS

STORAGE	Refer to the label
SHELF LIFE	2 Years

Analyst

Checker

Manager

Product Analysis Report


Product Name: 6"-Feruloylspinisin Chrom-Column: C18

Lot Number: 18061421 Mobile Phase: ACN:H2O=20:80

Detector: UV Wavelength: 254nm

Detect Time: 2018-06-14 Inject Vol: 20ul

Sample 50%MeOH Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	19.792	353.654	11159.222	0.1299
2	20.792	177.125	4284.979	0.0499
3	21.876	252745.734	8575432.000	99.8202