


PRODUCT ID	E-0717	Structure 
PRODUCT	Gomisin G/戈米辛G	
STANDARD TYPE	Qualitative Reference Material/含量测定	
LOT NUMBER	18030722	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-03-07	
DATE OF REPORT	2018-03-07	

OTHER NAME	NA
CHEMICAL FORMULA	C ₂₉ H ₃₂ O ₈
MOLECULAR WEIGHT(MW)	508.57
CAS NUMBER	62956-48-3
CHEMICAL FAMILY	Lignans

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK577—2008	NLT ≥98.0%	98.8%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK577—2008	NA	ND
Water	Q/TFFK577—2008	≤2%	Conforms
Appearance		White Crystal	Conforms

STORAGE CONDITIONS

STORAGE	Refer to the label
SHELF LIFE	2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Gomisin G Chrom-Column: C18
Lot Number: 18030722 Mobile Phase: MeOH:ACN:H2O=30:35:35
Detector: UV Wavelength: 254nm
Detect Time: 2018-03-07 Inject Vol: 20ul
Sample: MeOH Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	18.691	508333.250	15564142.000	98.8955
2	27.262	4033.164	173831.563	1.1045