



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-2444	<p align="center">Structure</p> 
PRODUCT	L-Dicentrine/ 荷包牡丹碱	
STANDARD TYPE	Qualitative Reference Material/含量测定	
LOT NUMBER	18010924	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-01-09	
DATE OF REPORT	2018-01-09	

OTHER NAME NA

CHEMICAL FORMULA C₂₀H₂₁NO₄

MOLECULAR WEIGHT(MW) 339.39

CAS NUMBER 28832-07-7

CHEMICAL FAMILY

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK337—2007	NLT ≥98.0%	99.9%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK337—2007	NA	ND
Water	Q/TFFK337—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE Refer to the label

SHELF LIFE 2 Years

Analyst

Checker

Manager

Product Analysis Report


Product Name: L-Dicentrine Chrom-Column: C18

Lot Number: 18010924 Mobile Phase: MeOH:0.1%H₃PO₄=40:60

Detector: UV Wavelength: 307nm

Detect Time: 2018-01-09 Inject Vol: 20ul

Sample 50%MeOH Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	6.235	407.298	8034.689	0.0837
2	10.237	239342.500	9590037.000	99.9163