



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

MSN:Tauto_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0703	Structure 
PRODUCT	Escin IA/七叶皂苷A	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	18051422	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-05-14	
DATE OF REPORT	2018-05-14	

OTHER NAME Aescin IA
CHEMICAL FORMULA C₅₅H₈₆O₂₄
MOLECULAR WEIGHT(MW) 1131.26
CAS NUMBER 123748-68-5
CHEMICAL FAMILY Saponin

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK671-2010	NLT ≥97.0%	97.1%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK671-2010	NA	ND
Water	Q/TFFK671-2010	≤2%	Conforms
Appearance		white powder	Conforms

STORAGE CONDITIONS

STORAGE Refer to the label
 SHELF LIFE 2 Years

Analyst

Checker

Manager

