



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-3110	<p style="text-align: center;">Structure</p> 
PRODUCT	REBAUDIOSIDE B/甜菊苷B	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	18032721	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-3-27	
DATE OF REPORT	2018-3-27	

OTHER NAME	NA
CHEMICAL FORMULA	C38H60O18
MOLECULAR WEIGHT(MW)	804.8722
CAS NUMBER	58543-17-2
CHEMICAL FAMILY	Diterpenoids

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK20—2004	NLT ≥98.0%	100.0%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK20—2004	NA	ND
Water	Q/TFFK20—2004	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE	Refer to the label
SHELF LIFE	2 Years

Analyst

Checker

Manager

Product Analysis Report

Product Name: REBAUDIOSIDE B **Chrom-Column:** Arcus C18


Lot Number: 18032721 **Mobile Phase:** ACN: H2O(35: 65)

Detect Time: 2018-3-27 **Detector:** ELSD

Inject Vol: 20.0 μ L **Tube Temp:** 110 $^{\circ}$ C

Flow: 1.0ml/min **Gas Flow:** 2.5L/min

Sample dissolution MeOH **Impactor:** Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	10.732	307537.656	3926120.750	100.0000